

产品名称: Cyclophilin B (5F10) Mouse Monoclonal Antibody
产品货号: AMM03622



产品概述 (Summary)

产品名称 (Production Name)	Cyclophilin B (5F10) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPIB
别名 (Alternative Names)	PPIB; CYPB; Peptidyl-prolyl cis-trans isomerase B; PPlase B; CYP-S1; Cyclophilin B; Rotamase B; S-cyclophilin; SCYLP
基因 ID (Gene ID)	5479
蛋白 ID (SwissProt ID)	P23284.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 24 kDa

研究背景 (Background)

产品名称: Cyclophilin B (5F10) Mouse Monoclonal Antibody
产品货号: AMM03622




PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human Hepatocarcinoma using Cyclophilin B (5F10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemical analysis of paraffin-embedded Human tonsils using Cyclophilin B (5F10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



产品名称: Cyclophilin B (5F10) Mouse Monoclonal Antibody
产品货号: AMM03622



Western blot analysis of Cyclophilin B (5F10) in PC-3, T47D, COS7, C6 lysates using Cyclophilin B (5F10) antibody



Western blot analysis of Cyclophilin B (5F10) in Hela, COS7, PC-3, Human placenta, Jurkat lysates using Cyclophilin B (5F10) antibody.

注意事项 (Note)

For research use only .