

产品名称: Transcription Initiation Factor IIB (2F6) Mouse
Monoclonal Antibody
产品货号: AMM03639



产品概述 (Summary)

产品名称 (Production Name)	Transcription Initiation Factor IIB (2F6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTF2B General transcription factor IIB; General transcription factor TFIIB; gtf2b; RNA
别名 (Alternative Names)	polymerase II transcription factor IIB; S300 II; S300-II; TF IIB; TF2B; TF2B_HUMAN; TFIIB; Transcription initiation factor IIB.
基因 ID (Gene ID)	2959
蛋白 ID (SwissProt ID)	Q00403.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa; Observed MW: 35 kDa

**产品名称: Transcription Initiation Factor IIB (2F6) Mouse
Monoclonal Antibody**
产品货号: AMM03639




研究背景 (Background)

General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of TFIIIB in Jurkat, MCF-7, COS7, 3T3 and Hela lysates using TFIIIB antibody.

注意事项 (Note)

For research use only .