

产品名称: CHD3 (7F6) Mouse Monoclonal Antibody
产品货号: AMM03643



产品概述 (Summary)

产品名称 (Production Name)	CHD3 (7F6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHD3 CHD3; Chromodomain-helicase-DNA-binding protein 3; CHD-3; ATP-dependent helicase CHD3; Mi-2 autoantigen 240 kDa protein; Mi2-alpha; Zinc finger helicase; hZFH
别名 (Alternative Names)	
基因 ID (Gene ID)	1107
蛋白 ID (SwissProt ID)	Q12873.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 227 kDa; Observed MW: 260 kDa

产品名称: CHD3 (7F6) Mouse Monoclonal Antibody
产品货号: AMM03643




研究背景 (Background)

Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Required for anchoring centrosomal pericentrin in both interphase and mitosis, for spindle organization and centrosome integrity.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of CHD3 (Cterminus) in NTERA-2, 3T3, C6, K562 and HeLa lysates using CHD3 (Cterminus) antibody.



Immunocytochemistry analysis of CHD3 (7F6) in HeLa using CHD3 (Cterminus) antibody.



产品名称: CHD3 (7F6) Mouse Monoclonal Antibody
产品货号: AMM03643



Immunoprecipitation analysis of CHD3 (7F6) in Hela lysates using CHD3 antibody.

注意事项 (Note)

For research use only .