

产品名称: JAK1 (8B8) Mouse Monoclonal Antibody
产品货号: AMM03662



产品概述 (Summary)

产品名称 (Production Name)	JAK1 (8B8) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	JAK1
别名 (Alternative Names)	JAK1; JAK1A; JAK1B; Tyrosine-protein kinase JAK1; Janus kinase 1; JAK-1
基因 ID (Gene ID)	3716
蛋白 ID (SwissProt ID)	P23458.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 133 kDa; Observed MW: 130 kDa

研究背景 (Background)

Tyrosine kinase of the non-receptor type, involved in the IFN-alpha/beta/gamma signal pathway. Kinase partner for the

产品名称: JAK1 (8B8) Mouse Monoclonal Antibody
产品货号: AMM03662




interleukin (IL)-2 receptor.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of JAK1 (8B8) in C6, Ramos, Jurkat and HeLa lysates using Jak1 antibody.



Immunocytochemistry analysis of JAK1 (8B8) in HeLa using Jak1 antibody.



Immunoprecipitation analysis of JAK1 (8B8) in HeLa lysates using Jak1 antibody.

产品名称: JAK1 (8B8) Mouse Monoclonal Antibody
产品货号: AMM03662



注意事项 (Note)

For research use only .