

产品名称: 14-3-3 eta (5F2) Mouse Monoclonal Antibody
产品货号: AMM03674



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 eta (5F2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAH 14-3-3 protein eta; 1433F_HUMAN; Brain protein 14-3-3; eta isoform;
别名 (Alternative Names)	HGNC:12853; Protein AS1; Tyrosine 3/tryptophan 5 monooxygenase activation protein eta polypeptide; YWHA 1; YWHA1; Ywhah.
基因 ID (Gene ID)	7533
蛋白 ID (SwissProt ID)	Q04917.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

产品名称: 14-3-3 eta (5F2) Mouse Monoclonal Antibody
产品货号: AMM03674




研究背景 (Background)

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of 14-3-3 eta/YWHAH in HeLa, 3T3, COS7 and C2C12 lysates using 14-3-3 eta/YWHAH antibody.

注意事项 (Note)

For research use only .