

产品名称: LIN28A (9G2) Mouse Monoclonal Antibody
产品货号: AMM03705



产品概述 (Summary)

产品名称 (Production Name)	LIN28A (9G2) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Ascitic Fluid

免疫原信息 (Immunogen)

基因名 (Gene Name)	LIN28A
别名 (Alternative Names)	CSDD1; FLJ12457; LIN 28; Lin-28A; LIN28; LIN28A; LN28A_HUMAN; Protein lin-28 homolog A; ZCCHC1; Zinc finger CCHC domain-containing protein 1.
基因 ID (Gene ID)	79727
蛋白 ID (SwissProt ID)	Q9H9Z2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

Involved in the timing of developmental events and choice of stage specific cell fates. Acts as a suppressor of microRNA

产品名称: LIN28A (9G2) Mouse Monoclonal Antibody
产品货号: AMM03705




(miRNA) biogenesis by specifically binding the precursor let-7 (pre-let-7), a miRNA precursor. Acts by binding pre-let-7 and recruiting ZCCHC11/TUT4 uridylyltransferase, leading to the terminal uridylation of pre-let-7.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Immunofluorescence analysis of LIN28A (9G2) in HeLa cells were transfected with pMX construct of Human LIN28 using LIN28A (9G2) antibody



Immunofluorescence analysis of LIN28A (9G2) in NTERA-2 cells using LIN28A (9G2) antibody (green). Blue: DRAQ5 fluorescent DNA dye.



产品名称: LIN28A (9G2) Mouse Monoclonal Antibody
产品货号: AMM03705



Western blot analysis of LIN28A (9G2) in NTERA-2 lysates using LIN28A (9G2) antibody.

注意事项 (Note)

For research use only .