

产品名称: GSK3 beta Mouse Monoclonal Antibody
产品货号: AMM03710



产品概述 (Summary)

产品名称 (Production Name)	GSK3 beta Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSK3B
别名 (Alternative Names)	GSK3B; Glycogen synthase kinase-3 beta; GSK-3 beta; Serine/threonine-protein kinase GSK3B
基因 ID (Gene ID)	2932
蛋白 ID (SwissProt ID)	P49841.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

产品名称: GSK3 beta Mouse Monoclonal Antibody
产品货号: AMM03710




Glycogen synthase kinase-3 (GSK3) is a proline-directed serine-threonine kinase that was initially identified as a phosphorylating and inactivating glycogen synthase. GSK3B is involved in energy metabolism, neuronal cell development, and body pattern formation. In skeletal muscle, it contributes to insulin regulation of glycogen synthesis by phosphorylating and inhibiting GYS1 activity and hence glycogen synthesis.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)




Flow Cytometry analysis of NIH/3T3 cells stained with GSK3B antibody (green) and negative control (red).




Immunohistochemistry analysis of paraffin-embedded rectum cancer tissues using GSK3B antibody with DAB staining. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: GSK3 beta Mouse Monoclonal Antibody
产品货号: AMM03710



Immunohistochemistry analysis of paraffin-embedded esophageal cancer tissues using GSK3B antibody with DAB staining. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of GSK3B antibody in HeLa, 3T3, C6, rat Brain lysates using GSK3B antibody.

注意事项 (Note)

For research use only .