

产品名称: GRP78 BiP (4F11) Mouse Monoclonal Antibody
产品货号: AMM03871



产品概述 (Summary)

产品名称 (Production Name)	GRP78 BiP (4F11) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPA5 HSPA5; GRP78; 78 kDa glucose-regulated protein; GRP-78; Endoplasmic reticulum lumenal Ca(2+)-binding protein grp78; Heat shock 70 kDa protein 5; Immunoglobulin heavy chain-binding protein; BiP
别名 (Alternative Names)	
基因 ID (Gene ID)	3309
蛋白 ID (SwissProt ID)	P11021.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 72 kDa; Observed MW: 78 kDa

产品名称: GRP78 BiP (4F11) Mouse Monoclonal Antibody
产品货号: AMM03871




研究背景 (Background)

When Chinese hamster K12 cells are starved of glucose, the synthesis of several proteins, called glucose-regulated proteins (GRPs), is markedly increased. Hendershot et al. (1994) (PubMed 8020977) pointed out that one of these, GRP78 (HSPA5), also referred to as 'immunoglobulin heavy chain-binding protein' (BiP), is a member of the heat-shock protein-70 (HSP70) family and is involved in the folding and assembly of proteins in the endoplasmic reticulum (ER).


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)




Western blot analysis of GRP78 BiP (4F11) in 3T3, A549, C6 and HeLa lysates using BiP/GRP78 antibody.



Immunofluorescence analysis of GRP78 BiP (4F11) in HeLa using BiP/GRP78 antibody.

产品名称: GRP78 BiP (4F11) Mouse Monoclonal Antibody
产品货号: AMM03871



Immunoprecipitation analysis of GRP78 BiP (4F11) in HeLa lysates using BiP/GRP78 antibody. Western blot analysis of BiP/GRP78 using BiP/GRP78 antibody.

注意事项 (Note)

For research use only .