

产品名称: EGFR (6H11) Mouse Monoclonal Antibody
产品货号: AMM03872



产品概述 (Summary)

产品名称 (Production Name)	EGFR (6H11) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EGFR
别名 (Alternative Names)	EGFR; ERBB; ERBB1; HER1; Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1
基因 ID (Gene ID)	1956
蛋白 ID (SwissProt ID)	P00533.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 134 kDa; Observed MW: 175 kDa

研究背景 (Background)

产品名称: EGFR (6H11) Mouse Monoclonal Antibody
产品货号: AMM03872




EGFR is a receptor tyrosine kinase. Receptor for epidermal growth factor (EGF) and related growth factors including TGF- α , amphiregulin, betacellulin, heparin-binding EGF-like growth factor, GP30 and vaccinia virus growth factor. Is involved in the control of cell growth and differentiation. . A single-pass transmembrane tyrosine kinase. Ligand binding to this receptor results in receptor dimerization, autophosphorylation (in trans), activation of various downstream signaling molecules and lysosomal degradation.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Western blot analysis of EGFR in A549 and COS7 lysates using EGFR antibody.




Immunocytochemistry analysis of EGFR (6H11) in HeLa using EGFR antibody.


产品名称: EGFR (6H11) Mouse Monoclonal Antibody
产品货号: AMM03872



Immunoprecipitation analysis of EGFR (6H11) in HeLa lysates using EGFR antibody.



Immunocytochemistry analysis of EGFR in MDA-MB-468 cells using EGFR antibody.



Western blot analysis of EGFR in A549 and MDA-MB-468 lysates using EGFR antibody.

注意事项 (Note)

For research use only .