

产品名称: Bcl2 (3H5) Mouse Monoclonal Antibody
产品货号: AMM03878



产品概述 (Summary)

产品名称 (Production Name)	Bcl2 (3H5) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Other,Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCL2
别名 (Alternative Names)	BCL2; Apoptosis regulator Bcl-2
基因 ID (Gene ID)	596
蛋白 ID (SwissProt ID)	P10415.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

研究背景 (Background)

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as

产品名称: Bcl2 (3H5) Mouse Monoclonal Antibody
产品货号: AMM03878



lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using Bcl2 (3H5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Western blot analysis of Bcl2 (3H5) in chicken lysates using Bcl2 (3H5) antibody



Western blot analysis of Bcl2 in Human breast cancer cell line MCF-7(A), MDAMB231(B) and Cal51(C) using Bcl2 antibody.

产品名称: Bcl2 (3H5) Mouse Monoclonal Antibody
产品货号: AMM03878



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Bcl2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .