

产品名称: Phosphoserine(5B12)Mouse Monoclonal Antibody
产品货号: AMM06003



产品概述 (Summary)

产品名称 (Production Name)	Phosphoserine(5B12)Mouse Monoclonal Antibody
描述 (Description)	Mouse Monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Species independent

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phospho Antibody
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, pH 7.4, containing 0.5%BSA, 0.02% New type preservative N as Preservative and 50% Glycerol.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	.Purified Protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:40000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

Serine (phosphatidylserine, PS), also known as serine phospholipids, two acyl serine. Phosphate ester compounds in the


产品名称: Phosphoserine(5B12)Mouse Monoclonal Antibody
产品货号: AMM06003




glycerol phosphate esters, cell membrane group, is one of the important constituents of cell membranes in the brain, especially in the nervous system of the body and also various functions of the brain (especially on the brain's memory and emotional stability) plays important roles, such as it can affects the membrane fluidity and permeability, and can activate the metabolism and synthesis of a variety of enzymes.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of 1) Hela, 2) Rat Brain Tissue, 3) Mouse Brain Tissue with Phosphoserine Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using Phosphoserine Mouse mAb diluted at 1:200.

产品名称: Phosphoserine(5B12)Mouse Monoclonal Antibody
产品货号: AMM06003



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Phosphoserine Mouse mAb diluted at 1:200.

注意事项 (Note)

For research use only .