

产品名称: Bovine Serum Albumin Mouse Monoclonal Antibody(Mix-mA)
产品货号: AMM07633



产品概述 (Summary)

产品名称 (Production Name)	Bovine Serum Albumin Mouse Monoclonal Antibody(Mix-mA)
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Cow

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, pH 7.4, containing 0.02% New type preservative N as Preservative and 50% Glycerol.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	.Purified Protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	66kDa

研究背景 (Background)


产品名称: Bovine Serum Albumin Mouse Monoclonal Antibody(Mix-mA)
产品货号: AMM07633



Albumin is the most abundant protein in the circulatory system. In mammals albumin is synthesized initially as preproalbumin by the liver. After removal of the signal peptide, the resultant proalbumin is further processed by removal of the six-residue propeptide from the new N-terminus. The albumin released into circulation possesses a half-life of 19 days.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Bovine Serum Albumin using Bovine Serum Albumin Monoclonal Antibody.

注意事项 (Note)

For research use only .