

产品名称: Chromogranin A (ABT037) Mouse Monoclonal Antibody
产品货号: AMM08783



产品概述 (Summary)

产品名称 (Production Name)	Chromogranin A (ABT037) Mouse Monoclonal Antibody
描述 (Description)	Mouse Monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ELISA
种属反应性 (Reactivity)	Human, Mouse, not react with Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHGA beta Granin;betagranin (N-terminal fragment of chromogranin A);catestatin;CgA;CHG A;Chga;chromofungin;Chromogranin A;Chromogranin A parathyroid secretory protein 1;Chromogranin A precursor;ChromograninA;CMGA_HUMAN;ER-37;Pancreastatin;Parastatin;Parathyroid secretory protein 1;Pituitary secretory protein I;Secretory protein I;SP I;SP-I;SP1;SPI;vasostatin 2;Vasostatin;Vasostatin I;Vasostatin II;vasostatin-2
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P10645.

产品应用 (Application)

产品名称: Chromogranin A (ABT037) Mouse Monoclonal Antibody
产品货号: AMM08783


稀释比 (Dilution Ratio)	IHC 1:200-400
蛋白分子量 (Molecular Weight)	50kD

研究背景 (Background)

function: Pancreastatin strongly inhibits glucose induced insulin release from the pancreas., miscellaneous: Binds calcium with a low-affinity., PTM: Sulfated on tyrosine residues and/or contains sulfated glycans., similarity: Belongs to the chromogranin/secretogranin protein family., subcellular location: Neuroendocrine and endocrine secretory granules.,

研究领域 (Research Area)

图片 (Image Data)



Human pancreas tissue was stained with Anti-Chromogranin A (ABT037) Antibody

注意事项 (Note)

For research use only .