

产品名称: Rabbit IgG (Fc specific, Heavy Chain)(1E1)Mouse
Mouse Monoclonal Antibody
产品货号: AMM16811



产品概述 (Summary)

产品名称 (Production Name)	Rabbit IgG (Fc specific, Heavy Chain)(1E1)Mouse Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Rabbit

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, pH 7.4, containing 0.5% protective protein, 0.02% New type preservative N as Preservative and 50% Glycerol.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	.Recombinant Protein of Rabbit IgG (Fc specific, Heavy Chain)

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:5000
蛋白分子量 (Molecular Weight)	55kDa

**产品名称: Rabbit IgG (Fc specific, Heavy Chain)(1E1)Mouse
Mouse Monoclonal Antibody
产品货号: AMM16811**




研究背景 (Background)

Whole IgG antibodies are isolated as intact molecules from antisera by immunoaffinity chromatography. They have an Fc portion and two antigen binding Fab portions joined together by disulfide bonds and therefore they are divalent. The average molecular weight is reported to be about 160 kDa. The Fc specific, Heavy Chain antibody reacts with the Fc portion of IgG heavy chain but not with the Fab portion of IgG.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Reduced Rabbit IgG, H+L, diluted at 1:10000.

注意事项 (Note)

For research use only .