

产品名称: CSF3 (8Q6) Mouse Monoclonal Antibody
产品货号: AMM20366



产品概述 (Summary)

产品名称 (Production Name)	CSF3 (8Q6) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS (pH 7.3) containing 1% protective protein, 50% glycerol and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CSF3
别名 (Alternative Names)	G-CSF; GCSF; CSF3OS; C17orf33
基因 ID (Gene ID)	1440.0
蛋白 ID (SwissProt ID)	P09919.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:200-1:500
蛋白分子量 (Molecular Weight)	21.4kDa

研究背景 (Background)

Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the

产品名称: CSF3 (8Q6) Mouse Monoclonal Antibody
产品货号: AMM20366




production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes macrophages. This CSF induces granulocytes.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Immunohistochemical staining of paraffin-embedded Human skin tissue within the normal limits using anti-CSF3 mouse monoclonal antibody.



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CSF3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CSF3 (1:2000).

注意事项 (Note)

For research use only .