

产品名称: ASB2 (5M19) Mouse Monoclonal Antibody
产品货号: AMM20414



产品概述 (Summary)

产品名称 (Production Name)	ASB2 (5M19) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS (pH 7.3) containing 1% protective protein, 50% glycerol and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity Purified

免疫原信息 (Immunogen)

基因名 (Gene Name)	ASB2
别名 (Alternative Names)	-
基因 ID (Gene ID)	NM_016150
蛋白 ID (SwissProt ID)	Q96Q27.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:20-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 64.9 kDa

研究背景 (Background)


产品名称: ASB2 (5M19) Mouse Monoclonal Antibody
产品货号: AMM20414




This gene encodes a member of the ankyrin repeat and SOCS box-containing (ASB) protein family. These proteins play a role in protein degradation by coupling suppressor of cytokine signalling (SOCS) proteins with the elongin BC complex. The encoded protein is a subunit of a multimeric E3 ubiquitin ligase complex that mediates the degradation of actin-binding proteins. This gene plays a role in retinoic acid-induced growth inhibition and differentiation of myeloid leukemia cells. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2011]

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemical staining of paraffin-embedded Carcinoma of Human thyroid tissue using anti-ASB2 mouse monoclonal antibody.




Immunohistochemical staining of paraffin-embedded Human thyroid tissue within the normal limits using anti-ASB2 mouse monoclonal antibody.

产品名称: ASB2 (5M19) Mouse Monoclonal Antibody
产品货号: AMM20414



Immunohistochemical staining of paraffin-embedded Carcinoma of Human prostate tissue using anti-ASB2 mouse monoclonal antibody.



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ASB2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ASB2 (1:500).

注意事项 (Note)

For research use only .