

产品名称: Iba 1 Mouse Monoclonal Antibody
产品货号: AMM20473



产品概述 (Summary)

产品名称 (Production Name)	Iba 1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IBA1
别名 (Alternative Names)	AIF 1 ;AIF-1 ;Aif1;BART 1 ;G1 ;IBA 1 ;IBA1 ;interferon gamma responsive transcript ;I ;IRT 1 ;MRF1 ;Protein g1 ;
基因 ID (Gene ID)	199.0
蛋白 ID (SwissProt ID)	P55008.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC/IF 1:100-1:500,ELISA 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 17 kDa

研究背景 (Background)

This gene encodes a protein that binds actin and calcium. This gene is induced by cytokines and interferon and may


产品名称: Iba 1 Mouse Monoclonal Antibody
产品货号: AMM20473




promote macrophage activation and growth of vascular smooth muscle cells and T-lymphocytes. Polymorphisms in this gene may be associated with systemic sclerosis. Alternative splicing results in multiple transcript variants, but the full-length and coding nature of some of these variants is not certain. [provided by RefSeq, Jan 2016],

研究领域 (Research Area)

图片 (Image Data)



Various whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Iba 1 antibody. The HRP-conjugated Goat anti-Mouse IgG(H + L) antibody was used to detect the antibody.



Human tonsil tissue was stained with Anti-Iba 1 Antibody

注意事项 (Note)

For research use only .