

产品名称: Apoa5 Mouse Monoclonal Antibody
产品货号: AMM80504



产品概述 (Summary)

产品名称 (Production Name)	Apoa5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Apoa5
别名 (Alternative Names)	RAP3; APOAV; APOA-V
基因 ID (Gene ID)	116519.0
蛋白 ID (SwissProt ID)	Q6Q788.Purified recombinant fragment of human APOA5 (AA: 20-363) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	41.2kDa

研究背景 (Background)

Apolipoprotein A5 (ApoA5) is fast gaining attention as a key regulator of serum triglyceride concentrations. An ApoA5

产品名称: Apoa5 Mouse Monoclonal Antibody
产品货号: AMM80504



mouse knock-out model produced an approximately four fold increase in serum triglycerides, whereas a knock-in model with human ApoA5 produced 50

研究领域 (Research Area)

图片 (Image Data)




Figure 1: Western blot analysis using Apoa5 mouse mAb against human serum (1) and Apoa5 recombinant protein (2).

注意事项 (Note)

For research use only .