

产品名称: CD4 Mouse Monoclonal Antibody
产品货号: AMM80505



产品概述 (Summary)

产品名称 (Production Name)	CD4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD4
别名 (Alternative Names)	CD4mut
基因 ID (Gene ID)	920.0
蛋白 ID (SwissProt ID)	P01730.Purified recombinant fragment of human CD4 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)

The CD4 antigen is highly expressed on human T helper cells and thymocytes, and at lower levels on monocytes and dendritic cells. This integral membrane glycoprotein of approximately 58 kDa contains four external domains (D1 to D4)


产品名称: CD4 Mouse Monoclonal Antibody
产品货号: AMM80505



that show homology to members of the immunoglobulin (Ig) superfamily. It is responsible for the recognition of the MHC class II antigen. The CD4 antibody recognizes most thymocytes and about 65% of peripheral blood T cells.

研究领域 (Research Area)

图片 (Image Data)



Flow cytometric analysis of blood T cells using CD4 mouse mAb (green) and negative control (blue).

注意事项 (Note)

For research use only .