

产品名称: KSHV K8a Mouse Monoclonal Antibody
产品货号: AMM80506



产品概述 (Summary)

产品名称 (Production Name)	KSHV K8a Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA
种属反应性 (Reactivity)	Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KSHV K8a
别名 (Alternative Names)	K8
基因 ID (Gene ID)	4961462.0
蛋白 ID (SwissProt ID)	.Purified recombinant fragment of KSHV K8α expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:25-1:50,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)


Kaposi's sarcoma-associated herpesvirus (KSHV) belongs to the gamma-(2)-herpesvirus subfamily and has been closely linked to the Kaposi

产品名称: KSHV K8a Mouse Monoclonal Antibody
产品货号: AMM80506




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using KSHV K8 mouse mAb against BCBL-1 (1) and TPA induced BCBL-1 (2) cell lysate.



Immunocytochemistry analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) using KSHV K8 mouse mAb with AEC staining.

注意事项 (Note)

For research use only .