

产品名称: KSHV ORF45 Mouse Monoclonal Antibody
产品货号: AMM80507



产品概述 (Summary)

产品名称 (Production Name)	KSHV ORF45 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	KSHV

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KSHV ORF45
别名 (Alternative Names)	ORF45
基因 ID (Gene ID)	4961474.0
蛋白 ID (SwissProt ID)	.Purified recombinant fragment of KSHV ORF45 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)


Kaposi's sarcoma-associated herpesvirus (KSHV) belongs to the gamma-(2)-herpesvirus subfamily and has been closely linked to the Kaposi

产品名称: KSHV ORF45 Mouse Monoclonal Antibody
产品货号: AMM80507



研究领域 (Research Area)

图片 (Image Data)



Immunocytochemistry analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) using KSHV ORF45 mouse mAb with AEC staining.

注意事项 (Note)

For research use only .