

产品名称: CD3 Mouse Monoclonal Antibody
产品货号: AMM80514



产品概述 (Summary)

产品名称 (Production Name)	CD3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD3
别名 (Alternative Names)	T3Z; CD247
基因 ID (Gene ID)	919.0
蛋白 ID (SwissProt ID)	P20963.Purified recombinant fragment of human CD3 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)

CD3 is a component of T cell surface receptor (TCR), it contains 5 invariable chains: gamma-, delta-, epsilon-, and zeta/eta. Clevers et al. (1988) showed the TCR/CD3 complex consists of either a TCR alpha/beta or TCR gamma/delta heterodimer.


产品名称: CD3 Mouse Monoclonal Antibody
产品货号: AMM80514



coexpressed at the cell surface with the invariant subunits,CD3 is responsible for the signal transduction of the TCR.

研究领域 (Research Area)

图片 (Image Data)



FCM analysis CD3 expressed in human peripheral blood cells.

注意事项 (Note)

For research use only .