

产品名称: BSA Mouse Monoclonal Antibody
产品货号: AMM80518



产品概述 (Summary)

产品名称 (Production Name)	BSA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BSA
别名 (Alternative Names)	BSA
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	.Bovine serum albumin.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 20,35 kDa

研究背景 (Background)

Bovine serum albumin (BSA) is used in variety of laboratories products which contain BSA as a
stabilizer. BSA could be detected for the indication of remains of chemical derivates from animal
products. This antibody is a good


产品名称: BSA Mouse Monoclonal Antibody
产品货号: AMM80518



reagent that may be used for the development of detection assay in ELISA or western blot to test the remaning BSA for clinical products in Vivo or research products In Vitro.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using BSA mouse mAb antiobdy against BSA.

注意事项 (Note)

For research use only .