

产品名称: APOA5 Mouse Monoclonal Antibody
产品货号: AMM80526



产品概述 (Summary)

产品名称 (Production Name)	APOA5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APOA5 (Pair)
别名 (Alternative Names)	RAP3; APOAV; APOA-V; MGC126836; MGC126838
基因 ID (Gene ID)	116519.0
蛋白 ID (SwissProt ID)	Q6Q788.Purified recombinant fragment of human Apoa5 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)

Apolipoprotein A5 (ApoA5) is fast gaining attention as a key regulator of serum triglyceride concentrations. An ApoA5 mouse knock-out model produced an approximately four fold increase in serum triglycerides, whereas a knock-in model


产品名称: APOA5 Mouse Monoclonal Antibody
产品货号: AMM80526



with human ApoA5 produced 50.

研究领域 (Research Area)

图片 (Image Data)



Sandwich ELISA using antibody pair to detect APOA5 protein.

注意事项 (Note)

For research use only .