

产品名称: ABL2 Mouse Monoclonal Antibody
产品货号: AMM80584



产品概述 (Summary)

产品名称 (Production Name)	ABL2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABL2
别名 (Alternative Names)	ARG; ABL; FLJ22224; FLJ31718; FLJ41441
基因 ID (Gene ID)	27.0
蛋白 ID (SwissProt ID)	P42684.Purified recombinant fragment of ABL2 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)

ABL2(ARG, Abelson-related gene) is a cytoplasmic tyrosine kinase which is closely related to but distinct from ABL1. The similarity of the proteins includes the tyrosine kinase domains and extends amino-terminal to include the SH2 and SH3

产品名称: ABL2 Mouse Monoclonal Antibody
产品货号: AMM80584



domains. ABL2 is expressed in both normal and tumor cells. The ABL2 gene product is expressed as two variants bearing different amino termini, both approximately 12-kb in length. c-Abl shows both cytoplasmic and nuclear localization, c-Abl is involved in two different chromosomal translocations present in human leukemias, which generate Bcr-Abl and TEL-Abl.

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence staining of methanol-fixed HeLa cells using ABL2 mouse mAb showing cytoplasm localization.

注意事项 (Note)

For research use only .