

产品名称: Thioredoxin (TRX) Mouse Monoclonal Antibody
产品货号: AMM80586



产品概述 (Summary)

产品名称 (Production Name)	Thioredoxin (TRX) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Tag

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Thioredoxin (TRX)
别名 (Alternative Names)	Thioredoxin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	.Purified recombinant fusion protein with Thioredoxin (TRX) tag.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:500-1:1000,ICC 1:100-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	12kDa

研究背景 (Background)

Thioredoxin (TRX) is a small ubiquitous protein (MW12kDa) which is exist in a wide variety of prokaryotic and eukaryotic


产品名称: Thioredoxin (TRX) Mouse Monoclonal Antibody
产品货号: AMM80586



cells. Trx contains a redox active disulfide/dithiol group within the conserved Cys-Gly-Pro-Cys active site. This antibody is suitable for detecting fusion proteins which encode a Trx-Tag by immunoblotting and immunoprecipitation. The Monoclonal Antibody can detect a little Trx-Tag fusion proteins with negligible cross-reactivity with bacterial, insect, or mammalian lysates.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using Trx mouse mAb against various fusion protein with Trx tag.

注意事项 (Note)

For research use only .