

产品名称: LPA Mouse Monoclonal Antibody
产品货号: AMM80703



产品概述 (Summary)

产品名称 (Production Name)	LPA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LPA
别名 (Alternative Names)	LP; AK38; APOA; LPA
基因 ID (Gene ID)	4018.0
蛋白 ID (SwissProt ID)	P08519 .Purified recombinant fragment of LPA (1823 - 2013) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)

LPA: lipoprotein, Lp(a). Apo(a) is the main constituent of lipoprotein(a) (Lp(a)). It has serine proteinase activity and is able of


产品名称: LPA Mouse Monoclonal Antibody
产品货号: AMM80703



autoproteolysis. Inhibits tissue-type plasminogen activator 1.Lp(a) may be a ligand for megalin/gp 330.Apo(a) is known to be proteolytically cleaved, leading to the formation of the so-called mini-lp(a).Apo(a) fragments accumulate in atherosclerotic lesions,where they may promote thrombogenesis.O-glycosylation may limit the extent of proteolytic fragmentation.

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human Liver tissues using LPA mouse mAb

注意事项 (Note)

For research use only .