

产品名称: Myostatin Mouse Monoclonal Antibody
产品货号: AMM80711



产品概述 (Summary)

产品名称 (Production Name)	Myostatin Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Myostatin
别名 (Alternative Names)	Myostatin; GDF8; MSTN
基因 ID (Gene ID)	2660.0
蛋白 ID (SwissProt ID)	O14793.Purified recombinant fragment of Myostatin expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:50-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	42.7kDa

研究背景 (Background)

The protein is a member of the bone morphogenetic protein (BMP) family and the TGF-beta superfamily. This group of proteins is characterized by a polybasic proteolytic processing site which is cleaved to produce a mature protein containing


产品名称: Myostatin Mouse Monoclonal Antibody
产品货号: AMM80711




seven conserved cysteine residues. The members of this family are regulators of cell growth and differentiation in both embryonic and adult tissues. This gene is thought to encode a secreted protein which negatively regulates skeletal muscle growth.

研究领域 (Research Area)


图片 (Image Data)



Flow cytometric analysis of COS7 cells using Myostatin mouse mAb (green) and negative control (red).



Flow cytometric analysis of HeLa cells using Myostatin mouse mAb (green) and negative control (red).



Flow cytometric analysis of NIH3T3 cells using Myostatin mouse mAb (green) and negative control (red).

产品名称: Myostatin Mouse Monoclonal Antibody
产品货号: AMM80711



注意事项 (Note)

For research use only .