

产品名称: DKK1 Mouse Monoclonal Antibody
产品货号: AMM80720



产品概述 (Summary)

产品名称 (Production Name)	DKK1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DKK1
别名 (Alternative Names)	SK; DKK-1;
基因 ID (Gene ID)	22943.0
蛋白 ID (SwissProt ID)	O94907.Purified recombinant fragment of DKK1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	28.7kDa

研究背景 (Background)

DKK1: dickkopf homolog 1 (*Xenopus laevis*), also known as SK. Entrez Protein NP_036374. DKK1 is a member of the dickkopf family. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through

产品名称: DKK1 Mouse Monoclonal Antibody
产品货号: AMM80720




its inhibition of the WNT signaling pathway. Elevated levels of DKK1 in bone marrow plasma and peripheral blood is associated with the presence of osteolytic bone lesions in patients with multiple myeloma.


研究领域 (Research Area)

Wnt signaling pathway


图片 (Image Data)



Flow cytometric analysis of COS7 cells using DKK1 mouse mAb (green) and negative control (red).



Flow cytometric analysis of hela cells using DKK1 mouse mAb (green) and negative control (red).



Flow cytometric analysis of NIH3T3 cells using DKK1 mouse mAb (green) and negative control (red).

产品名称: DKK1 Mouse Monoclonal Antibody
产品货号: AMM80720



注意事项 (Note)

For research use only .