

产品名称: human IgG(Fc Specific) Mouse Monoclonal Antibody
产品货号: AMM80721



产品概述 (Summary)

产品名称 (Production Name)	human IgG(Fc Specific) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	human IgG(Fc Specific)
别名 (Alternative Names)	human IgG(Fc Specific)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	.Purified fusion protein with human IgG(Fc Specific) tag.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)

Specificity of antibody has been tested using Western blot, immunofluorescence, and ELISA. The product reacts with


产品名称: human IgG(Fc Specific) Mouse Monoclonal Antibody
产品货号: AMM80721




human IgG (Fc) with no significant reactivity with human IgG (Fab2), IgM or other serum proteins. We didn't test crossing with other species IgG (Fc).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using human IgG (Fc specific) mouse mAb against different fusion proteins with human IgG(Fc specific) tag.



Confocal Immunofluorescence analysis of HEK293 cells transfected with recombinant plasmid with human IgG Fc tag using human IgG Fc mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.

注意事项 (Note)

For research use only .