

产品名称: MPS1 Mouse Monoclonal Antibody
产品货号: AMM80745



产品概述 (Summary)

产品名称 (Production Name)	MPS1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MPS1
别名 (Alternative Names)	RPS27
基因 ID (Gene ID)	6232.0
蛋白 ID (SwissProt ID)	P42677.Purified recombinant fragment of MPS1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	9.5kDa

研究背景 (Background)

MPS1, also known as RPS27. It is a ribosomal protein. Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80


产品名称: MPS1 Mouse Monoclonal Antibody
产品货号: AMM80745



structurally distinct proteins. MPS1 is a component of the 40S subunit. The protein belongs to the S27E family of ribosomal proteins. It contains a C4-type zinc finger domain that can bind to zinc. The encoded protein has been shown to be able to bind to nucleic acid. It is located in the cytoplasm as a ribosomal component, but it has also been detected in the nucleus. Studies in rat indicate that ribosomal protein S27 is located near ribosomal protein S18 in the 40S subunit and is covalently linked to translation initiation factor eIF3. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

研究领域 (Research Area)

图片 (Image Data)



Confocal Immunofluorescence analysis of Hela cells (A), BCBL-1 cells (B) and L1210 cells (C) using MPS1 mouse mAb (green).
Blue: DRAQ5 fluorescent DNA dye.

注意事项 (Note)

For research use only .