

产品名称: TBX5 Mouse Monoclonal Antibody
产品货号: AMM80748



产品概述 (Summary)

产品名称 (Production Name)	TBX5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBX5
别名 (Alternative Names)	HOS; TBX5
基因 ID (Gene ID)	6910.0
蛋白 ID (SwissProt ID)	Q99593.Purified recombinant fragment of TBX5 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	53kDa

研究背景 (Background)

TBX5, also known as T-box 5. It is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental


产品名称: TBX5 Mouse Monoclonal Antibody
产品货号: AMM80748



processes. It is closely linked to related family member T-box 3 (ulnar mammary syndrome) on human chromosome 12. The TBX5 protein may play a role in heart development and specification of limb identity. Mutations in this gene have been associated with Holt-Oram syndrome, a developmental disorder affecting the heart and upper limbs. Several transcript variants encoding different isoforms have been described for this gene.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using TBX5 mouse mAb against HepG2 cell lysate (1).

注意事项 (Note)

For research use only .