

产品名称: PDGFR Mouse Monoclonal Antibody
产品货号: AMM80751



产品概述 (Summary)

产品名称 (Production Name)	PDGFR Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDGFR
别名 (Alternative Names)	JTK12; PDGFR; CD140B; PDGFR1; PDGFRB
基因 ID (Gene ID)	5159.0
蛋白 ID (SwissProt ID)	P09619.Purified recombinant fragment of human PDGFR β expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	190kDa

研究背景 (Background)

Platelet-derived growth factor receptor, beta polypeptide, also known as PDGFRB, PDGFR


产品名称: PDGFR Mouse Monoclonal Antibody
产品货号: AMM80751



研究领域 (Research Area)

PI3K-Akt signaling pathway, MAPK signaling pathway, Hippo signaling pathway

图片 (Image Data)



Western blot analysis using PDGFR mouse mAb against NIH/3T3 cell lysate (1).

注意事项 (Note)

For research use only .