

产品名称: Neurogenin-3 Mouse Monoclonal Antibody
产品货号: AMM80767



产品概述 (Summary)

产品名称 (Production Name)	Neurogenin-3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Neurogenin-3
别名 (Alternative Names)	ngn3; Atoh5; Math4B; bHLHa7; NEUROG3
基因 ID (Gene ID)	50674.0
蛋白 ID (SwissProt ID)	Q9Y4Z2.Purified recombinant fragment of human Neurogenin-3 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	23kDa

研究背景 (Background)

Neurogenin-3 (NEUROG3), also known as NGN3, it is expressed in endocrine progenitor cells and is required for endocrine


产品名称: Neurogenin-3 Mouse Monoclonal Antibody
产品货号: AMM80767




cell development in the pancreas and intestine. It belongs to a family of basic helix-loop-helix transcription factors involved in the determination of neural precursor cells in the neuroectoderm.

研究领域 (Research Area)


图片 (Image Data)



Flow cytometric analysis of HeLa cells using Neurogenin-3 mouse mAb (green) and negative control (red).



Flow cytometric analysis of NIH3T3 cells using Neurogenin-3 mouse mAb (green) and negative control (red).



Flow cytometric analysis of PC-12 cells using Neurogenin-3 mouse mAb (green) and negative control (red).

产品名称: Neurogenin-3 Mouse Monoclonal Antibody
产品货号: AMM80767



注意事项 (Note)

For research use only .