

产品名称: TEC Mouse Monoclonal Antibody
产品货号: AMM80768



产品概述 (Summary)

产品名称 (Production Name)	TEC Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TEC
别名 (Alternative Names)	PSCTK4
基因 ID (Gene ID)	7006.0
蛋白 ID (SwissProt ID)	P42680.Purified recombinant fragment of TEC expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	74kDa

研究背景 (Background)

TEC, it belongs to the Tec family of non-receptor protein-tyrosine kinases containing a pleckstrin homology domain. Tec family kinases are involved in the intracellular signaling mechanisms of cytokine receptors, lymphocyte surface antigens,


产品名称: TEC Mouse Monoclonal Antibody
产品货号: AMM80768




heterotrimeric G-protein coupled receptors, and integrin molecules. They are also key players in the regulation of the immune functions. Tec kinase is an integral component of T cell signaling and has a distinct role in T cell activation. This gene may be associated with myelodysplastic syndrome.

研究领域 (Research Area)

图片 (Image Data)



Flow cytometric analysis of Jurkat cells using TEC mouse mAb (green) and negative control (red).



Flow cytometric analysis of NIH3T3 cells using TEC mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .