

产品名称: FRK Mouse Monoclonal Antibody
产品货号: AMM80769



产品概述 (Summary)

产品名称 (Production Name)	FRK Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FRK
别名 (Alternative Names)	GTK; RAK; PTK5
基因 ID (Gene ID)	2444.0
蛋白 ID (SwissProt ID)	P42685 .Purified recombinant fragment of human FRK expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	58kDa

研究背景 (Background)

FRK (Fyn-related kinase), also known as Rak. It is an epithelial tissue-specific kinase. The human FRK gene maps to chromosome 6q21-q22.3 and encodes a 505 amino acid protein. The protein belongs to the TYR family of protein kinases.


产品名称: FRK Mouse Monoclonal Antibody
产品货号: AMM80769



This tyrosine kinase is a nuclear protein and may function during G1 and S phase of the cell cycle and suppress growth.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using FRK mouse mAb against K562 cell lysate (1).

注意事项 (Note)

For research use only .