

产品名称: GTF2B Mouse Monoclonal Antibody
产品货号: AMM80793



产品概述 (Summary)

产品名称 (Production Name)	GTF2B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTF2B
别名 (Alternative Names)	TF2B; TFIIB; GTF2B
基因 ID (Gene ID)	2959.0
蛋白 ID (SwissProt ID)	Q00403.Purified recombinant fragment of human GTF2B expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	34.8kDa

研究背景 (Background)

This gene encodes the general transcription factor IIB, one of the ubiquitous factors required for transcription initiation by RNA polymerase II. The protein localizes to the nucleus where it forms a complex (the DAB complex) with transcription


产品名称: GTF2B Mouse Monoclonal Antibody
产品货号: AMM80793



factors IID and IIA. Transcription factor IIB serves as a bridge between IID, the factor which initially recognizes the promoter sequence, and RNA polymerase II.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using GTF2B mouse mAb against HeLa (1), NIH/3T3 (2), COS7 (3) and A431 (4) cell lysate.

注意事项 (Note)

For research use only .