

产品名称: PRMT6 Mouse Monoclonal Antibody
产品货号: AMM80798



产品概述 (Summary)

产品名称 (Production Name)	PRMT6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRMT6
别名 (Alternative Names)	HRMT1L6; FLJ10559; FLJ51477; PRMT6
基因 ID (Gene ID)	55170.0
蛋白 ID (SwissProt ID)	Q96LA8.Purified recombinant fragment of human PRMT6 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	42kDa

研究背景 (Background)

Protein arginine N-methyltransferases, such as PRMT6, catalyze the sequential transfer of a methyl group from S-adenosyl-L-methionine to the side chain nitrogens of arginine residues within proteins to form methylated arginine derivatives and


产品名称: PRMT6 Mouse Monoclonal Antibody
产品货号: AMM80798



S-adenosyl-L-homocysteine.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using PRMT6 mouse mAb against A431 (1), HeLa (2), A549 (3) and HEK293 (4) cell lysate.

注意事项 (Note)

For research use only .