

产品名称: TBLR1 Mouse Monoclonal Antibody
产品货号: AMM80807



产品概述 (Summary)

产品名称 (Production Name)	TBLR1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBLR1
别名 (Alternative Names)	C21; DC42; IRA1; TBL1XR1
基因 ID (Gene ID)	79718.0
蛋白 ID (SwissProt ID)	Q9BZK7.Purified recombinant fragment of human TBLR1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	55.5kDa

研究背景 (Background)

The protein encoded by this gene has sequence similarity with members of the WD40 repeat-containing protein family. The WD40 group is a large family of proteins, which appear to have a regulatory function. It is believed that the WD40 repeats


产品名称: TBLR1 Mouse Monoclonal Antibody
产品货号: AMM80807




mediate protein-protein interactions and members of the family are involved in signal transduction, RNA processing, gene regulation, vesicular trafficking, cytoskeletal assembly and may play a role in the control of cytotypic differentiation.

研究领域 (Research Area)


图片 (Image Data)



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Western blot analysis using TBLR1 mouse mAb against Mouse heart (1) tissue lysate and K562 (2) cell lysate.




Flow cytometric analysis of K562 cells using TBLR1 mouse mAb (green) and negative control (purple).

产品名称: TBLR1 Mouse Monoclonal Antibody
产品货号: AMM80807



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using TBLR1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using TBLR1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .