

产品名称: CSF1 Mouse Monoclonal Antibody
产品货号: AMM80816



产品概述 (Summary)

产品名称 (Production Name)	CSF1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CSF1
别名 (Alternative Names)	MCSF
基因 ID (Gene ID)	1435.0
蛋白 ID (SwissProt ID)	P09603.Purified recombinant fragment of human CSF1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)

CSF1, also known as MCSF, is a four- α -helical bundle cytokine that controls the production, differentiation, and function of macrophages. The active form of the protein is found extracellularly as a disulfide-linked homodimer, and is thought to be

产品名称: CSF1 Mouse Monoclonal Antibody
产品货号: AMM80816




produced by proteolytic cleavage of membrane-bound precursors. This protein may be involved in development of the placenta.

研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human Placenta tissues using anti-CSF1 mouse mAb

注意事项 (Note)

For research use only .