

产品名称: ETS1 Mouse Monoclonal Antibody
产品货号: AMM80826



产品概述 (Summary)

产品名称 (Production Name)	ETS1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ETS1
别名 (Alternative Names)	ETS-1; EWSR2; FLJ10768; ETS1
基因 ID (Gene ID)	2113.0
蛋白 ID (SwissProt ID)	P14921.Purified recombinant fragment of human ETS1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	50.4kDa

研究背景 (Background)

The ETS1 is a transcription factor known to regulate expression of a number of genes involved in extracellular matrix


产品名称: ETS1 Mouse Monoclonal Antibody
产品货号: AMM80826




remodeling. ETS1 is absent from normal gastric epithelium but is expressed in approximately 60% of gastric carcinomas and oral squamous cell carcinomas with a significant correlation to the tumor stage.

研究领域 (Research Area)


图片 (Image Data)



Immunofluorescence analysis of NIH/3T3 cells using anti-ETS1 mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of Jurkat cells using anti-ETS1 mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human Colon, Adipocytes tissues using anti-ETS1 mouse mAb

产品名称: ETS1 Mouse Monoclonal Antibody
产品货号: AMM80826



注意事项 (Note)

For research use only .