

产品名称: PRDM1 Mouse Monoclonal Antibody
产品货号: AMM80828



产品概述 (Summary)

产品名称 (Production Name)	PRDM1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRDM1
别名 (Alternative Names)	BLIMP1; PRDI-BF1; MGC118922; MGC118923; MGC118924; MGC118925; PRDM1
基因 ID (Gene ID)	639.0
蛋白 ID (SwissProt ID)	O75626.Purified recombinant fragment of human PRDM1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	88kDa

研究背景 (Background)

PRDM1/Blimp1 is a repressor of beta-interferon gene expression. This action is exerted by binding to the PRDI (positive


产品名称: PRDM1 Mouse Monoclonal Antibody
产品货号: AMM80828



regulatory domain I element) of the beta-IFN gene promoter. The transcription of this gene is increased upon virus induction. Two alternatively spliced transcript variants that encode different isoforms have been reported.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using PRDM1 mouse mAb against Raji (1, 2), L1210 (3) and TPH-1 (4) cell lysate.

注意事项 (Note)

For research use only .