

产品名称: MUM1 Mouse Monoclonal Antibody
产品货号: AMM80848



产品概述 (Summary)

产品名称 (Production Name)	MUM1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MUM1
别名 (Alternative Names)	MUM-1; EXPAND1; HSPC211; FLJ14868; FLJ22283; MGC131891; MGC163315; MUM1
基因 ID (Gene ID)	84939.0
蛋白 ID (SwissProt ID)	Q2TAK8.Purified recombinant fragment of human MUM1 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	78.6kDa

研究背景 (Background)

Mutated melanoma-associated antigen 1 (MUM1) is a 78 kDa protein encoded by the MUM1 gene. MUM1 is expressed in


产品名称: MUM1 Mouse Monoclonal Antibody
产品货号: AMM80848




a subset of B cells in the light zone of the germinal center, plasma cells, activated T cells and a wide spectrum of related hematolymphoid neoplasms. Antibodies to MUM1 are useful in the sub classification of lymphoid malignancies.

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human melanoma (A), non-Hodgkin's lymphoma (B), using MUM1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human Uterus tissues using MUM1 mouse mAb

注意事项 (Note)

For research use only .