

产品名称: KLHL12 Mouse Monoclonal Antibody
产品货号: AMM80855



产品概述 (Summary)

产品名称 (Production Name)	KLHL12 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLHL12
别名 (Alternative Names)	DKIR; C3IP1; FLJ27152; KLHL12
基因 ID (Gene ID)	59349.0
蛋白 ID (SwissProt ID)	Q53G59.Purified recombinant fragment of human KLHL12 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	63kDa

研究背景 (Background)

Serves as a substrate-specific adapter for the CUL3-based ubiquitin-protein E3 ligase complex. Negatively regulates the Wnt signaling pathway via the targeted ubiquitination and subsequent proteolysis of DVL3. Tissue specificity: Highly

产品名称: KLHL12 Mouse Monoclonal Antibody
产品货号: AMM80855




expressed in testis and at lower levels in the submandibular salivary gland.

研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Western blot analysis using KLHL12 mouse mAb against Hela (1) cell lysate.

注意事项 (Note)

For research use only .