

产品名称: KLHL21 Mouse Monoclonal Antibody
产品货号: AMM80857



产品概述 (Summary)

产品名称 (Production Name)	KLHL21 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLHL21
别名 (Alternative Names)	KIAA0469; MGC99635; KLHL21
基因 ID (Gene ID)	9903.0
蛋白 ID (SwissProt ID)	Q9UJP4.Purified recombinant fragment of human KLHL21 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	66kDa

研究背景 (Background)


Probable substrate-specific adapter of an E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

产品名称: KLHL21 Mouse Monoclonal Antibody
产品货号: AMM80857



研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human brain tissues using KLHL21 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .