

产品名称: KLHL11 Mouse Monoclonal Antibody
产品货号: AMM80862



产品概述 (Summary)

产品名称 (Production Name)	KLHL11 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLHL11
别名 (Alternative Names)	FLJ10572; KLHL11
基因 ID (Gene ID)	55175.0
蛋白 ID (SwissProt ID)	Q9NVR0 .Purified recombinant fragment of human KLHL11 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	80kDa

研究背景 (Background)


Kelch-like protein 11, Subcellular location Secreted.

产品名称: KLHL11 Mouse Monoclonal Antibody
产品货号: AMM80862




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using KLHL11 mouse mAb against Hela (1) and MCF-7 (2) cell lysate.



Immunofluorescence analysis of 3T3-L1 cells using KLHL11 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of K562 cells using KLHL11 mouse mAb (green) and negative control (purple).

产品名称: KLHL11 Mouse Monoclonal Antibody
产品货号: AMM80862



注意事项 (Note)

For research use only .