

产品名称: IL34 Mouse Monoclonal Antibody
产品货号: AMM80864



产品概述 (Summary)

产品名称 (Production Name)	IL34 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL34
别名 (Alternative Names)	IL-34; C16orf77; MGC34647; IL34
基因 ID (Gene ID)	146433.0
蛋白 ID (SwissProt ID)	Q6ZMJ4 .Purified recombinant fragment of human IL34 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	27.5kDa

研究背景 (Background)

Cytokine that promotes the differentiation and viability of monocytes and macrophages. Stimulates phosphorylation of MAPK1/ERK2 AND MAPK3/ERK1. Ligand for colony-stimulating factor-1 receptor CSF1R.Tissue specificity: Detected in the


产品名称: IL34 Mouse Monoclonal Antibody
产品货号: AMM80864




sinusoidal epithelium in the red pulp of spleen (at protein level). Predominantly expressed in spleen. Also detected in a range of other tissues including heart, brain, lung, liver, kidney, thymus, testis, ovary, small intestine, prostate and colon.

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of U251 cells using IL34 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of K562 cells using IL34 mouse mAb (green) and negative control (purple).

注意事项 (Note)

For research use only .