

产品名称: KLF4 Mouse Monoclonal Antibody
产品货号: AMM80874



产品概述 (Summary)

产品名称 (Production Name)	KLF4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLF4
别名 (Alternative Names)	EZF; GKLF; KLF4
基因 ID (Gene ID)	9314.0
蛋白 ID (SwissProt ID)	O43474.Purified recombinant fragment of human KLF4 (AA: 2-470) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

研究背景 (Background)

Kruppel-like factor 4 (gut),endothelial Kruppel-like zinc finger protein. May act as a transcriptional activator. Binds the


产品名称: KLF4 Mouse Monoclonal Antibody
产品货号: AMM80874



CACCC core sequence. May be involved in the differentiation of epithelial cells and may also function in the development of the skeleton and kidney. KLF4 is highly expressed in the epithelial cells of the skin and the gastrointestinal tract.

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues (left) and human rectum tissues (right) using KLF4 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .