

产品名称: NFKBIB Mouse Monoclonal Antibody
产品货号: AMM80875



产品概述 (Summary)

产品名称 (Production Name)	NFKBIB Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFKBIB
别名 (Alternative Names)	IKBB; TRIP9; NFKBIB
基因 ID (Gene ID)	4793.0
蛋白 ID (SwissProt ID)	Q15653.Purified recombinant fragment of human NFKBIB expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)

Inhibits NF-kappa-B by complexing with and trapping it in the cytoplasm. However, the unphosphorylated form resynthesized after cell stimulation is able to bind NF-kappa-B allowing its transport to the nucleus and protecting it to


产品名称: NFKBIB Mouse Monoclonal Antibody
产品货号: AMM80875



further IKBA-dependent inactivation. Association with inhibitor kappa B-interacting NKIRAS1 and NKIRAS2 prevent its phosphorylation rendering it more resistant to degradation, explaining its slower degradation. Tissue specificity: Expressed in all tissues examined.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using NFKBIB mouse mAb against Jurkat (1) cell lysate.

注意事项 (Note)

For research use only .