

产品名称: THAP11 Mouse Monoclonal Antibody
产品货号: AMM80876



产品概述 (Summary)

产品名称 (Production Name)	THAP11 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 0.03% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THAP11
别名 (Alternative Names)	RONIN; CTG-B43a; CTG-B45d; HRIHFB2206; THAP11
基因 ID (Gene ID)	57215.0
蛋白 ID (SwissProt ID)	Q96EK4.Purified recombinant fragment of human THAP11 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	44kDa

研究背景 (Background)

The protein encoded by this gene contains a THAP domain, which is a conserved DNA-binding domain that has striking


产品名称: THAP11 Mouse Monoclonal Antibody
产品货号: AMM80876




similarity to the site-specific DNA-binding domain (DBD) of Drosophila P element transposases.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using THAP11 mouse mAb against Hela (1) and NTERA-2 (2) cell lysate.




Immunofluorescence analysis of NTERA-2 cells using THAP11 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of Hela cells using THAP11 mouse mAb (green) and negative control (purple).

产品名称: THAP11 Mouse Monoclonal Antibody
产品货号: AMM80876



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues (left) and ovary cancer tissues (right) using THAP11 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .